



THE JAS EZ SHOULDER ADVANCED TURNBUCKLE TECHNOLOGY

The JAS EZ Shoulder offers full-range, multi-plane Range of Motion (ROM) therapy – 100° external rotation, 90° internal rotation and 90° abduction – in a lightweight, low-profile, single-patient use device. An innovative scapular stabilization feature, combined with infinitely adjustable ROM, assures precise capsular stretch throughout the entire ROM. Maximal Total End Range Time (TERT) is achieved for optimal therapeutic benefit.



JAS QUALITY AND PERFORMANCE IN A SINGLE-PATIENT USE DESIGN

**Full-Range,
Multi-Plane
Shoulder ROM
Therapy**

**NEW!
Body Cuff
Option Now
Available**

**Optimal
Fit and
Compliance**

**Patented
Motion
Arm™**

The **JAS EZ Shoulder** is complete with a fully adjustable **floor stand** or **body cuff**, custom cuffing and a scapular stabilization component that assures simple fit and accurate shoulder alignment every treatment session.

JAS EZ Shoulder

100° External Rotation
90° Internal Rotation
90° Abduction

Effective for Use Following:

- Adhesive Capsulitis
- Rotator Cuff Tear/Repairs
- Humeral Fractures
- Bankhart Procedures
- Mastectomy
- Total/Hemi Shoulder Replacements

AVAILABLE FOR COMMERCIAL PURCHASE



Joint Active Systems, Inc. / Telephone: 217.342.3412 or 800.879.0117 / Info@jointactivesystems.com


www.jointactivesystems.com

Licensed in the State of Illinois.



SHOULDER

SUPERIOR DESIGN. PROVEN OUTCOMES.



**Scapular
Stabilization
Feature**



**NEW!
Body Cuff Option
Allows Patient
Mobility**



**Easy-Fit
Floor Stand
Option
Quick and
Simple to
Apply**



**Custom
Cuffing for
Optimal Fit and
Performance**

Scapular Stabilization Feature

A unique scapular stabilization component maximizes stretch precisely to capsular tissues of the gleno-humeral joint.

Comfortable Body Cuff Option

A custom-fit body cuff offers versatility, portability, and comfort, and assures accurate shoulder alignment during use, maximizing compliance and outcomes.

Convenient Floor Stand Option

A fully adjustable floor stand design makes it simple and comfortable for patients to apply, and assures accurate shoulder alignment during use, maximizing compliance and outcomes.

Custom Cuffing

Custom cuffing assures an accurate fit with no slipping or gapping, and maximizes soft tissue capture to optimize delivery of stretch loads to the tissues.

**Full-Range, Multi-Plane ROM Therapy –
100° External Rotation,
90° Internal Rotation and 90° Abduction**

AVAILABLE FOR COMMERCIAL PURCHASE



Joint Active Systems, Inc. / Telephone: 217.342.3412 or 800.879.0117 / Info@jointactivesystems.com

www.jointactivesystems.com

Licensed in the State of Illinois.